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Name	Type	Lot ID	Exp. date
APTT-SP	Intermediate reagent	N0917362	09/2013
APTT-SP CaCl2	Start reagent	N0917375	09/2014
AT Liq FXa	Intermediate reagent		
AT Liq Substrate	Start reagent		
Clean B	Clean		
Clean B Diluted	Clean		
D-Dimer High Control	Quality Control	B11789	02/2014
D-Dimer H5 Buffer	Intermediate reagent	B11780	04/2013
D-Dimer H5 Cal	Calibrator/NPP	B11780	04/2013
D-Dimer H5 Latex	Start reagent	B11780	04/2013
D-Dimer Low Control	Quality Control	B11789	02/2014
Factor Diluent	Sample Diluent		
Fib-C	Start reagent		
Fib-C Abn Control	Quality Control		
FXIII Buffer	Intermediate reagent		
FXIII Diluent	Sample Diluent		
FXIII Latex	Start reagent		
HemosIL Cal Plasma	Calibrator/NPP	E0616328	06/2014
High Abn C Unassayed	Quality Control		
High Abn C. Assayed	Quality Control	N0616330	06/2014
Low Abn C. Assayed	Quality Control	N1118082	11/2014
Low Abn C. Unassayed	Quality Control		
Low Fib Control	Quality Control		
Normal C. Assayed	Quality Control	N1118092	11/2014
Normal C. Unassayed	Quality Control		

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General Information		Stir, Rinse & Clean Information		Lot Specific Information	
Material index:	681	Category:	Reagents		
Material name:	APTT TS	Material type:	Intermediate Reagent		
Manufacturer:	TS	Bottle type:	10mL		
Volume tracking <input checked="" type="checkbox"/> Enable volume threshold: 0.4 mL		On-board stability <input checked="" type="checkbox"/> Enable on-board stability tracking On-board stability: 48 Hours <input checked="" type="checkbox"/> Enable warning threshold: 1 Hours <input checked="" type="checkbox"/> Test not feasible when stability is expired		Comments:	
<input checked="" type="checkbox"/> Enable expiration tracking <input checked="" type="checkbox"/> Enable warning threshold: 30 Days					
Location information					
Track Identifier		Rack Identifier		Rack Pos.Identifier	
				Lot ID	

General Information		Stir, Rinse & Clean Information		Lot Specific Information	
Material index:	661	Category:	Reagents		
Material name:	APTT TS CaCl2	Material type:	Start Reagent		
Manufacturer:	TS	Bottle type:	10mL		
Volume tracking <input checked="" type="checkbox"/> Enable volume threshold: 0.4 mL		On-board stability <input checked="" type="checkbox"/> Enable on-board stability tracking On-board stability: 48 Hours <input checked="" type="checkbox"/> Enable warning threshold: 1 Hours <input checked="" type="checkbox"/> Test not feasible when stability is expired		Comments:	
<input checked="" type="checkbox"/> Enable expiration tracking <input checked="" type="checkbox"/> Enable warning threshold: 30 Days					
Location information					
Track Identifier		Rack Identifier		Rack Pos.Identifier	
				Lot ID	

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LIS number: 132 Test code: PT-RP(E)

Load cycle
Number of replicates: 1

Material Name	Type	Load Volume	Incubation Range
Sample		50	55,000 - 360,000
PT RecombiPlasTin 2G	Start Reagent	100	

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Material/Sample Load Cycle Definition

Material: **Sample**

Aspiration information

Head volume: 0 μL

Airgap: 15 μL

Head volume airgap: 0 μL

Material/Sample volume: 50 μL

Transport airgap: 30 μL

Total volume: 95 μL

Dispensation information

☐ Enhanced dispense

☐ Enable mix

Mix: 50 %

Rinse time after mix: 1 s

Number of cycles: 1

☒ Enable incubation time

Incubation range: **60.000 - 60.000 s**

Rinse & clean information

☒ Rinse

Time: 1 s

☐ Enable agitation: 1

☐ Clean & rinse

Clean & rinse cycles: 1

Clean material: System clean

Aspiration cycles: 1

Clean airgap: 15 μL

Clean volume: 130 μL

Clean transport airgap: 0 μL

Clean total volume: 145 μL

Hold time: 0

☐ Enable agitation: 1

Clean volume (all cycles): 130 μL

Frequency: Between changes in material only

Rinse time: 1 s

☐ Enable agitation: 1

☐ Lock clean & rinse

OK Cancel

Material/Sample Load Cycle Definition

Material: **APTT TS**

Aspiration information

Head volume: 0 μL

Airgap: 15 μL

Head volume airgap: 0 μL

Material/Sample volume: 50 μL

Transport airgap: 30 μL

Total volume: 95 μL

Dispensation information

☐ Enhanced dispense

☐ Enable mix

Mix: 50 %

Rinse time after mix: 1 s

Number of cycles: 1

☒ Enable incubation time

Incubation range: **180.000 - 180.000 s**

Rinse & clean information

☒ Rinse

Time: 1 s

☐ Enable agitation: 1

☐ Clean & rinse

Clean & rinse cycles: 1

Clean material: System clean

Aspiration cycles: 1

Clean airgap: 15 μL

Clean volume: 130 μL

Clean transport airgap: 0 μL

Clean total volume: 145 μL

Hold time: 0

☐ Enable agitation: 1

Clean volume (all cycles): 130 μL

Frequency: Between changes in material only

Rinse time: 1 s

☐ Enable agitation: 1

☐ Lock clean & rinse

OK Cancel

Material/Sample Load Cycle Definition

Material: **CACL2 TS**

Aspiration information

Head volume: 0 µL

Airgap: 15 µL

Head volume airgap: 0 µL

Material/Sample volume: **50** µL

Transport airgap: 30 µL

Total volume: 95 µL

Dispensation information

☐ Enhanced dispense

☐ Enable mix

Mix: 50 %

Rinse time after mix: 1 s

Number of cycles: 1

☒ Enable incubation time

Incubation range: 20,000 - 150,000 s

Rinse & clean information

☒ Rinse

Time: 1 s

☐ Enable agitation: 1

☐ Clean & rinse

Clean & rinse cycles: 1

Clean material: System clean

Aspiration cycles: 1

Clean airgap: 15 µL

Clean volume: 130 µL

Clean transport airgap: 0 µL

Clean total volume: 145 µL

Hold time: 0

☐ Enable agitation: 1

Clean volume (all cycles): 130 µL

Frequency: Between changes in material only

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Test Definition - APTT TS

- General information
- Analytical cycle definition
 - Sample pre-dilution
 - Alternative pre-dilution
- DR Parameters
 - Primary wavelength
 - Wavelength definition
 - Raw data checks
 - Normalized data checks
 - Primary algorithm
 - Secondary algorithm
- Result units**
- Calibration
 - Automatic dilution
 - Pre-diluted calibrators
- Parallelism
- Normal Pool Plasma
- Rerun - Rules

Result unit definition

Unit	Unit Type	Label	dec.	Normal Range	Therapeutic Range	Linear Range	Test Range
Seconds	Measured	s	<input checked="" type="checkbox"/>	1	25.0 36.5		16.0 500.0
Ratio	Calculated	Ratio	<input checked="" type="checkbox"/>	2	0.00 1.20		
			<input type="checkbox"/>				
			<input type="checkbox"/>				
			<input type="checkbox"/>				
			<input type="checkbox"/>				
			<input type="checkbox"/>				

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
ACL TOP Family - Version 4.4.0 - [Material List]								
Instrument Actions Analysis QC Calibration NPP Setup System Help								
Name Type Lot ID Exp. date								
APTT-SP	Intermediate reagent	N0917362	09/2013					
APTT-SP CaCl2	Start reagent	N0917375	09/2014					
AT Liq FXa	Intermediate reagent							
AT Liq Substrate	Start reagent							
Clean B	Clean							
Clean B Diluted	Clean							
D-Dimer High Control	Quality Control	B11789	02/2014					
D-Dimer HS Buffer	Intermediate reagent	B11780	04/2013					
D-Dimer HS Cal	Calibrator/NPP	B11780	04/2013					
D-Dimer HS Latex	Start reagent	B11780	04/2013					
D-Dimer Low Control	Quality Control	B11789	02/2014					
Factor Diluent	Sample Diluent							
Fib-C	Start reagent							
Fib-C Abn Control	Quality Control							
FXIII Buffer	Intermediate reagent							
FXIII Diluent	Sample Diluent							
FXIII Latex	Start reagent							
HemosIL Cal Plasma	Calibrator/NPP	E0616328	06/2014					
High Abn C Unassayed	Quality Control							
High Abn C. Assayed	Quality Control	N0616330	06/2014					
Low Abn C. Assayed	Quality Control	N1118082	11/2014					
Low Abn C. Unassayed	Quality Control							
Low Fib Control	Quality Control							
Normal C. Assayed	Quality Control	N1118092	11/2014					
Normal C. Unassayed	Quality Control							

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Instrument Actions Analysis QC Calibration NPP Setup System Help	
<div> </div>	
<div>Test Definition - FIB-C</div> <ul style="list-style-type: none"> General information Analytical cycle definition <ul style="list-style-type: none"> Sample pre-dilution Alternative pre-dilution DR Parameters <ul style="list-style-type: none"> Primary wavelength <ul style="list-style-type: none"> Wavelength definition Raw data checks Normalized data checks Primary algorithm Secondary algorithm Result units Calibration <ul style="list-style-type: none"> Automatic dilution Pre-diluted calibrators Parallelism Normal Pool Plasma Rerun - Rules 	<div>General information</div> <div> Test code: <input type="text" value="FIB-C"/> Test name: <input type="text" value="Fibrinogen Clauss"/> </div> <div> LIS number: <input type="text" value="1310"/> Test number: <input type="text" value="13"/> <input checked="" type="checkbox"/> IL test <input type="checkbox"/> For investigation only </div> <div>Assay</div> <div> Assay release number: <input type="text" value="7.4"/> </div> <div> Test modification number: <input type="text" value="220"/> </div> <div> User revision comment: <input type="text"/> </div> <div> IL revision comment: <input type="text"/> </div> <div> <input type="checkbox"/> Enable as shadow test </div> <div> Parent test: <input type="text"/> </div> <div> <input type="checkbox"/> Enable as paired test </div> <div> Master paired test: <input type="text"/> </div> <div> <input type="checkbox"/> Enable paired test dependency </div>

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Test Definition - FIB-C

General information

Test code: FIB-C Test name: Fibrinogen Clauss

LIS number: 1310 Test number: 13 ☒ IL test: ☐ For investigation only

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Load cycle

Number of replicates: 1

Material Name	Type	Load Volume	Incubation Range
Sample		50	55,000 - 360,000
FIB-C	Start Reagent	100	

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
Unit 1: Seconds
Unit 2: g/L
Unit 3: mg/dL
Unit 4:


Primary unit: Seconds
Max difference: 10 %

Measured result curve display settings
☒ Enable auto scale
☐ Raw
☒ Normalized (Absorbance)
 Y-Axis minimum: 0.000 Y-Axis maximum: 3000.000 mAbs

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Test code: FIB-C Test name: Fibrinogen Clauss

LIS number: 1310 Test number: 13 ☒ IL test ☐ For investigation only

Assay release number: 7.4
Test modification number: 220
User revision comment:
IL revision comment:

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Test Definition -

- General information
- Analytical cycle definition
 - Sample pre-dilution
 - Alternative pre-dilution
- DR Parameters
 - Primary wavelength
 - Wavelength definition
 - Raw data checks
 - Normalized data checks
 - Primary algorithm
 - Secondary algorithm
- Result units
- Calibration**
 - Automatic dilution
 - Pre-diluted calibrators
- Parallelism
- Normal Pool Plasma
- Rerun - Rules

☒ Enable calibration

General | Frequency | DR Parameters

General

Number of replicates: 4

Calibration unit: ng/mL

Calibration Mode

☐ Automatic dilutions process

☒ Pre-diluted calibrators

☐ Enable import calibration

Imported test:

Imported type: Raw Data

Correction Factor: 1.000

Imported unit:

Test Feasibility

☐ Test not feasible when no validated calibration

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ACL TOP Family - Version 4.4.0 - [Material List]							
Instrument Actions Analysis QC Calibration NPP Setup System Help							
	Name	Type					
	APTT-SP	Intermediate reagent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	APTT-SP CaCl2	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	N0917362
	AT Liq FXa	Intermediate reagent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	N0917375
	AT Liq Substrate	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	Clean B	Clean	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
	Clean B Diluted	Clean	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
	D-Dimer High Control	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B11789
	D-Dimer HS Buffer	Intermediate reagent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	B11780
	D-Dimer HS Cal	Calibrator/NPP	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	B11780
	D-Dimer HS Latex	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	B11780
	D-Dimer Low Control	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B11789
	Factor Diluent	Sample Diluent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
	Fib-C	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	Fib-C Abn Control	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	FXIII Buffer	Intermediate reagent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	FXIII Diluent	Sample Diluent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	FXIII Latex	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	HemosIL Cal Plasma	Calibrator/NPP	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	E0616328
	High Abn C Unassayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	High Abn C. Assayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	N0616330
	Low Abn C. Assayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	N1118082
	Low Abn C. Unassayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	Low Fib Control	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
	Normal C. Assayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	N1118092
	Normal C. Unassayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	

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General Information		Stir, Rinse & Clean Information	Lot Specific Information
Material index:	501	Category:	Reagents
Material name:	PT TS	Material type:	Start Reagent
Manufacturer:	TS	Bottle type:	10mL
Volume tracking <input checked="" type="checkbox"/> Enable volume threshold: 0.4 mL		On-board stability <input checked="" type="checkbox"/> Enable on-board stability tracking On-board stability: 48 Hours <input checked="" type="checkbox"/> Enable warning threshold: 1 Hours	
<input checked="" type="checkbox"/> Enable expiration tracking <input checked="" type="checkbox"/> Enable warning threshold: 30 Days		<input checked="" type="checkbox"/> Test not feasible when stability is expired	
Comments:			
Location Information			
Track Identifier	Rack Identifier	Rack Pos.Identifier	Lot ID

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ACL TOP Family - Version 4.4.0 - [IL Material Definition - Start reagent PT RecombiPlasTin]

Instrument Actions Analysis QC Calibration NPP Setup System Help

General Information Stir, Rinse & Clean Information **Lot Specific Information**

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Instrument Actions Analysis QC Calibration NPP **Setup** System Help

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Instrument Actions Analysis QC Calibration NPP Setup System Help	
<div> <div> Test Definition - PT-RP <div> General information Analytical cycle definition Sample pre-dilution Alternative pre-dilution DR Parameters Primary wavelength Wavelength definition Raw data checks Normalized data checks Primary algorithm Secondary algorithm Result units Calibration Automatic dilution Pre-diluted calibrators Parallelism Normal Pool Plasma Rerun - Rules </div> </div> </div>	
<div> <div>General information</div> <div> Test code: PT-RP Test name: PT RecombiPlasTin LIS number: 1310 Test number: 13 <input checked="" type="checkbox"/> IL test <input type="checkbox"/> For investigation only </div> <div> <div>Assay</div> <div> Assay release number: 7,4 Test modification number: 220 User revision comment: IL revision comment: </div> </div> <div> <input type="checkbox"/> Enable as shadow test <div> Parent test: </div> </div> <div> <input type="checkbox"/> Enable as paired test <div> Master paired test: </div> <input type="checkbox"/> Enable paired test dependency </div> </div>	

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Test Definition - PT-RP

- General information**
- Analytical cycle definition
- Sample pre-dilution
- Alternative pre-dilution
- DR Parameters
 - Primary wavelength
 - Wavelength definition
 - Raw data checks
 - Normalized data checks
 - Primary algorithm
 - Secondary algorithm
- Result units
- Calibration
 - Automatic dilution
 - Pre-diluted calibrators
- Parallelism
- Normal Pool Plasma
- Rerun - Rules

Test code: PT-RP Test name: PT RecombiPlasTin

LIS number: 1310 Test number: 13 ☒ IL test ☐ For investigation only

Assay
Assay release number: 7.4
Test modification number: 220
User revision comment:
IL revision comment:

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ACL TOP Family - Version 4.4.0 - [Material List]									
Instrument Actions Analysis QC Calibration NPP Setup System Help									
	Name	Type						Lot ID	Exp. date
	APTT-SP	Intermediate reagent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>		N0917362	09/2013
	APTT-SP CaCl2	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>		N0917375	09/2014
	AT Liq FXa	Intermediate reagent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	AT Liq Substrate	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	Clean B	Clean	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			
	Clean B Diluted	Clean	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			
	D-Dimer High Control	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B11789		02/2014
	D-Dimer HS Buffer	Intermediate reagent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	B11780		04/2013
	D-Dimer HS Cal	Calibrator/NPP	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	B11780		04/2013
	D-Dimer HS Latex	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	B11780		04/2013
	D-Dimer Low Control	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B11789		02/2014
	Factor Diluent	Sample Diluent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			
	Fib-C	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	Fib-C Abn Control	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	FXIII Buffer	Intermediate reagent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	FXIII Diluent	Sample Diluent	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	FXIII Latex	Start reagent	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	HemosIL Cal Plasma	Calibrator/NPP	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	E0616328		06/2014
	High Abn C Unassayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	High Abn C. Assayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	N0616330		06/2014
	Low Abn C. Assayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	N1118082		11/2014
	Low Abn C. Unassayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	Low Fib Control	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			
	Normal C. Assayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	N1118092		11/2014
	Normal C. Unassayed	Quality Control	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>			

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General Information		Stir, Rinse & Clean Information		Lot Specific Information									
Material index:	661	Category:	Reagents										
Material name:	Thrombin TS	Material type:	Start Reagent										
Manufacturer:	TS	Bottle type:	10mL										
Volume tracking <input checked="" type="checkbox"/> Enable volume threshold: 0.4 mL		On-board stability <input checked="" type="checkbox"/> Enable on-board stability tracking On-board stability: 48 Hours <input checked="" type="checkbox"/> Enable warning threshold: 1 Hours		Comments:									
<input checked="" type="checkbox"/> Enable expiration tracking <input checked="" type="checkbox"/> Enable warning threshold: 30 Days		<input checked="" type="checkbox"/> Test not feasible when stability is expired											
Location information <table border="1"> <thead> <tr> <th>Track Identifier</th> <th>Rack Identifier</th> <th>Rack Pos.Identifier</th> <th>Lot ID</th> </tr> </thead> <tbody> <tr> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>						Track Identifier	Rack Identifier	Rack Pos.Identifier	Lot ID				
Track Identifier	Rack Identifier	Rack Pos.Identifier	Lot ID										

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

Test Definition - TT5 General information Analytical cycle definition Sample pre-dilution Alternative pre-dilution DR Parameters Primary wavelength Wavelength definition Raw data checks Normalized data checks Primary algorithm Secondary algorithm Result units Calibration Automatic dilution Pre-diluted calibrators Parallelism Normal Pool Plasma Rerun - Rules		Acquisition - wavelength Primary wavelength: W-671 nm Acquisition - time Standard time: 100 s Delay time: 5 s <input checked="" type="checkbox"/> Enable extended test mode Extended time: 320 s LIS number: 132 Test code: Load cycle Number of replicates: 1 <table border="1"> <thead> <tr> <th>Material Name</th> <th>Type</th> <th>Load Volume</th> <th>Incubation Range</th> </tr> </thead> <tbody> <tr> <td>Sample</td> <td></td> <td>80</td> <td>55,000 - 360,000</td> </tr> <tr> <td>Thrombin</td> <td>Start Reagent</td> <td>80</td> <td></td> </tr> <tr> <td> </td> <td></td> <td></td> <td></td> </tr> <tr> <td> </td> <td></td> <td></td> <td></td> </tr> <tr> <td> </td> <td></td> <td></td> <td></td> </tr> <tr> <td> </td> <td></td> <td></td> <td></td> </tr> </tbody> </table> Total cuvette volume: 150 uL Consistency Check. (Double click to go to the inconsistency).	Material Name	Type	Load Volume	Incubation Range	Sample		80	55,000 - 360,000	Thrombin	Start Reagent	80																	
Material Name	Type	Load Volume	Incubation Range																											
Sample		80	55,000 - 360,000																											
Thrombin	Start Reagent	80																												

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Test Definition - TT TS

- General information
- Analytical cycle definition**
 - Sample pre-dilution
 - Alternative pre-dilution
- DR Parameters
 - Primary wavelength
 - Wavelength definition
 - Raw data checks
 - Normalized data checks
 - Primary algorithm
 - Secondary algorithm
- Result units
- Calibration
 - Automatic dilution
 - Pre-diluted calibrators
- Parallelism
- Normal Pool Plasma
- Rerun - Rules

Acquisition - wavelength
Primary wavelength: W-671 nm

Acquisition - time
Standard time: 100 s
Delay time: 5 s

☒ Enable extended test mode

Extended time: 320 s

LIS number: 132 Test code:

Load cycle
Number of replicates: 1

Material Name	Type	Load Volume	Incubation Range
Sample		80	55,000 - 360,000
Thrombin	Start Reagent	80	

Total cuvette volume: 150 uL

Consistency Check. (Double click to go to the inconsistency).

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Normal Pool Plasma

Result unit definition

Unit	Unit Type	Label	dec.	Normal Range	Therapeutic Range	Linear Range	Test Range
Seconds	Measured	s	<input checked="" type="checkbox"/>	1	10.3 16.6		5.0 300.0
Ratio	Calculated	Ratio	<input checked="" type="checkbox"/>	2			

Patient result selection

Unit 1: Seconds

Unit 2:

Unit 3:

Unit 4:

Primary unit cfg.

Primary unit: Seconds

☒ Enable max difference

Max difference: 10 %

Measured result curve display settings

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General Information		Stir, Rinse & Clean Information		Lot Specific Information	
Material index:	501	Category:	Reagents		
Material name:	AT SUBS TS	Material type:	Start Reagent		
Manufacturer:	TS	Bottle type:	10mL		
Volume tracking <input checked="" type="checkbox"/> Enable volume threshold: 0.4 mL		On-board stability <input checked="" type="checkbox"/> Enable on-board stability tracking On-board stability: 48 Hours <input checked="" type="checkbox"/> Enable warning threshold: 1 Hours <input checked="" type="checkbox"/> Test not feasible when stability is expired		Comments:	
<input checked="" type="checkbox"/> Enable expiration tracking <input checked="" type="checkbox"/> Enable warning threshold: 30 Days					
Location information					
Track Identifier		Rack Identifier		Rack Pos.Identifier	
Lot ID					

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

Instrument Actions Analysis QC Calibration NPP Setup System Help	
<div> </div>	
Test Definition - AT LIQ <ul style="list-style-type: none"> General information Analytical cycle definition <ul style="list-style-type: none"> Sample pre-dilution Alternative pre-dilution DR Parameters <ul style="list-style-type: none"> Primary wavelength <ul style="list-style-type: none"> Wavelength definition Raw data checks Normalized data checks Primary algorithm Secondary algorithm Result units Calibration <ul style="list-style-type: none"> Automatic dilution Pre-diluted calibrators Parallelism Normal Pool Plasma Rerun - Rules 	General information Test code: <input type="text" value="AT LIQ"/> Test name: <input type="text" value="Antithrombin Liquid"/> LIS number: <input type="text" value="2071"/> Test number: <input type="text" value="13"/> <input checked="" type="checkbox"/> IL test <input type="checkbox"/> For investigation only Assay Assay release number: <input type="text" value="7.4"/> Test modification number: <input type="text" value="220"/> User revision comment: <input type="text"/> IL revision comment: <input type="text"/> <input type="checkbox"/> Enable as shadow test Parent test: <input type="text"/> <input type="checkbox"/> Enable as paired test Master paired test: <input type="text"/> <input type="checkbox"/> Enable paired test dependency

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Test Definition - AT LIQ

- General information
- Analytical cycle definition**
 - Sample pre-dilution
 - Alternative pre-dilution
- DR Parameters
 - Primary wavelength
 - Wavelength definition
 - Raw data checks
 - Normalized data checks
 - Primary algorithm
 - Secondary algorithm
- Result units
- Calibration
 - Automatic dilution
 - Pre-diluted calibrators
 - Parallelism
 - Normal Pool Plasma
 - Rerun - Rules

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LIS number: 132 Test code: PT-RP(E)

Load cycle
Number of replicates: 1

Material Name	Type	Load Volume	Incubation Range
Sample		50	55,000 - 360,000
AT Liq Fxa	Int. Reagent	50	
AT Liq Substrate	Start Reagent	50	

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Test Definition - AT TS

- General information
- Analytical cycle definition
 - Sample pre-dilution**
 - Alternative pre-dilution
- DR Parameters
 - Primary wavelength
 - Wavelength definition
 - Raw data checks
 - Normalized data checks
 - Primary algorithm
 - Secondary algorithm
- Result units
- Calibration
 - Automatic dilution
 - Pre-diluted calibrators
- Parallelism
- Normal Pool Plasma
- Rein - Rules

☒ Enable sample pre-dilution

General | Sample/Mixture | Diluent

Dilution definition

Number of parts sample:

Number of parts diluent:

Minimum aspiratable volume: μL

Volumes (μL)		
Sample	Mixture	Diluent
4		360

☒ Enable mix

Mix volume: % Number of cycles:

Rinse after mix: s

Test Definition - AT TS

- General information
- Analytical cycle definition
 - Sample pre-dilution**
 - Alternative pre-dilution
- DR Parameters
 - Primary wavelength
 - Wavelength definition
 - Raw data checks
 - Normalized data checks
 - Primary algorithm
 - Secondary algorithm
- Result units
- Calibration
 - Automatic dilution
 - Pre-diluted calibrators
- Parallelism

☒ Enable sample pre-dilution

General | **Sample/Mixture** | Diluent

Volumes

Head volume: μL

Airgap: μL

Transport airgap: μL

Head volume airgap: μL

Rinse

Time: s

Test Definition - AT TS

- General information
- Analytical cycle definition
 - Sample pre-dilution**
 - Alternative pre-dilution
- DR Parameters
 - Primary wavelength
 - Wavelength definition
 - Raw data checks
 - Normalized data checks
 - Primary algorithm
 - Secondary algorithm
- Result units

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